

Technical Bulletin

Bulletin Number: 239

oxdot Revision Number: $oxdot{R-01}$

Release Date: September 9, 1981

Revision Date: July 30, 1992

TOMORROW'S MANUFACTURING LEADERSHIP

TODAY

SUBJECT:

OSCILLATOR

MODEL:

DNCB 630/360-5 Select and 6 Select

REASON:

increase the diameter of the pin the

oscillator increasing its strength.

BEGINNING SERIAL NUMBER:

0001-2470CF and up have the new

oscillator factory installed. (Approximate Date: 09/18/81)

Other parts will change at the same time, therefore a kit will be necessary when ordering replacement oscillators for venders serialized prior to 0001-2470CF.

KIT WILL CONSIST OF:	1	D801,200,990.41	Oscillator-3D
	1	A801,201,030.01	Linkage Arm
	1	801,804,770.01	Nyliner (small) For Linkage Arm
	2	901,700,650.31	Can Retainer
	1	801,803,170.01	Nyliner (large)
			For Oscillator

Oscillator Assembly

Can Retainer

Nyliner

Pin

View Of Oscillator (For Reference Only)